|  |            |                                  |
|--|------------|----------------------------------|
| Drew Scientific Ltd                        | <b>CE</b>  | Control Kit Lot Number M-0172076 |
| Glycated Haemoglobin Two Level Control Kit |            | Level 1 Lot Number MO-00172722   |
| Part Number 014-175                        | <b>IVD</b> | Level 2 Lot Number MO-00172723   |
|  |            | Expiry Date 2008AUG31            |

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**For in vitro diagnostic use only.** These controls are intended for use in a quality control procedure for determining imprecision and accuracy of glycated haemoglobin measurements.

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The routine use of Control Materials is recommended for the long term monitoring of imprecision of analytical methods and is a requirement of Good Laboratory Practice.

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### **Contents, Contenu du kit, Inhalt, Contenuto, Conteúdo, Contenido**

The kit contains three vials of each of two levels. Each vial contains a preparation of freeze dried whole blood.

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**Store unopened** at 2 – 8 °C. The expiration date is given at the top of this page. After reconstitution, the material is stable for 7 days at 2 – 8 °C.

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**Remove the stopper** carefully. Using a volumetric pipette, measure 0.3mL of deionised water and add to the vial to reconstitute the control. Replace the stopper firmly and mix gently. Allow to stand for 5 minutes. Visually inspect that all the material has dissolved before use. Treat the reconstituted control as whole blood for the purpose of preparing samples for your glycated haemoglobin analysis method.

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**Do not** use time-expired material. Do not use incompletely dissolved material.

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**All the** individual blood donor lots used to prepare this material were tested by FDA approved methods and found to be non-reactive for hepatitis B surface antigen (HbsAG), hepatitis C (HCV), HIV-1 and HIV-2. No test methods can offer complete assurance that these materials are free from infective agents. Handle, use and dispose of these control materials as if they were human specimens and potentially infective.

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### **Assigned values**

| <b>Methods</b>  | <b>Level 1 Value</b> | <b>Level 1 Guidance Range</b> | <b>Level 2 Value</b> | <b>Level 2 Guidance Range</b> |
|---|----------------------|-------------------------------|----------------------|-------------------------------|
| Drew Scientific <b>DS5</b><br>Hb A1c <b>NGSP Aligned</b>      | <b>5.3</b>           | <b>4.2 – 6.4</b>              | <b>10.6</b>          | <b>8.8 – 12.4</b>             |
| Drew Scientific <b>DS5</b><br>Hb A1c <b>IFCC Aligned</b>      | <b>3.4</b>           | <b>2.3 – 4.5</b>              | <b>9.3</b>           | <b>7.4 – 11.2</b>             |
| Drew Scientific <b>G15</b><br>Hb A1c <b>NGSP Aligned</b>      | <b>5.5</b>           | <b>4.8 – 6.2</b>              | <b>11.2</b>          | <b>10.1 – 12.3</b>            |
| Drew Scientific <b>G15</b><br>Hb A1c <b>IFCC Aligned</b>      | <b>3.6</b>           | <b>2.8 – 4.4</b>              | <b>9.9</b>           | <b>8.7 – 11.1</b>             |
| Drew Scientific <b>HbGold</b><br>HbA1c ( <b>Non Aligned</b> ) | <b>4.5</b>           | <b>3.7 – 5.3</b>              | <b>10.1</b>          | <b>8.7 – 11.5</b>             |