

**EX-TROL™**  
Hematology Control

CE IVD



**LOT** EX018



2008-04-05

Quality Control Months: Jan, Feb, Mar.

DREW H1800,18MS,18MS+	<b>LOT</b> EX018-Low				<b>LOT</b> EX018-Normal				<b>LOT</b> EX018-High			
	Mean	Limits	Range		Mean	Limits	Range		Mean	Limits	Range	
WBC K/uL	2.0	± 0.4	1.6	2.4	7.4	± 0.8	6.6	8.2	18.6	± 2.2	16.4	20.8
Lymp K/uL	1.1	± 0.3	0.8	1.4	2.1	± 0.5	1.6	2.6	2.2	± 1.1	1.1	3.3
Mid K/uL	0.3	± 0.1	0.2	0.4	0.5	± 0.3	0.2	0.8	0.7	± 0.6	0.1	1.3
Gran K/uL	0.7	± 0.2	0.5	0.9	4.7	± 0.6	4.1	5.3	15.7	± 1.5	14.2	17.2
Lymp %	54.9	± 7.0	47.9	61.9	29.0	± 5.0	24.0	34.0	11.7	± 5.0	6.7	16.7
Mid %	12.5	± 4.0	8.5	16.5	6.6	± 4.0	2.6	10.6	3.6	± 3.0	0.6	6.6
Gran %	32.5	± 6.0	26.5	38.5	64.4	± 6.0	58.4	70.4	84.6	± 6.0	78.6	90.6
RBC M/uL	2.32	± 0.15	2.17	2.47	4.64	± 0.20	4.44	4.84	5.57	± 0.24	5.33	5.81
HGB g/dL	6.0	± 0.4	5.6	6.4	13.7	± 0.6	13.1	14.3	17.9	± 0.8	17.1	18.7
HCT %	17.1	± 2.0	15.1	19.1	39.5	± 2.4	37.1	41.9	51.7	± 2.8	48.9	54.5
MCV fL	73.6	± 4.0	69.6	77.6	85.1	± 4.0	81.1	89.1	92.8	± 4.0	88.8	96.8
MCH pg	25.9	± 2.4	23.5	28.3	29.5	± 2.8	26.7	32.3	32.2	± 2.8	29.4	35.0
MCHC g/dL	35.2	± 3.0	32.2	38.2	34.6	± 3.0	31.6	37.6	34.7	± 3.0	31.7	37.7
RDW %	15.9	± 4.0	11.9	19.9	17.2	± 4.0	13.2	21.2	16.4	± 5.0	11.4	21.4
PLT K/uL	75	± 20	55	95	238	± 35	203	273	441	± 60	381	501
MPV fL	8.3	± 2.0	6.3	10.3	6.9	± 2.0	4.9	8.9	6.5	± 2.0	4.5	8.5

Platelet Interference alarms (or PLT/RBC overlap) and WBC differential "R" alarms may be seen.