



DREW Scientific, Inc.

Product No. 610020

Lot No. CDC010

Expiration : 5-Feb-2010

### MD CAL-KIT™

## HEMATOLOGY CALIBRATOR HEMAVET SERIES

For use with HEMAVET® HEMATOLOGY INSTRUMENTS

Parameter (units)	Target	Range
WBC (K/ $\mu$ L)	9.2	8.8 – 9.6
RBC (M/ $\mu$ L)	4.48	4.31 – 4.65
Hb (g/dL)	13.6	13.2 – 14.0
MCV (fL)	84	83.0 – 86.0
PLT (K/ $\mu$ L)	298	273 – 323

### IMPORTANT

Follow mixing instructions from reverse side and  
Wipe threads on bottle and cap before closing

**FOR IN VITRO DIAGNOSTIC USE  
PROTECT FROM FREEZING  
STORE AT 2° - 8° C (35° - 46° F)**

### POTENTIAL BIOHAZARDOUS MATERIAL

4230 Shilling Way Dallas, Texas. 75237 Phone (214) 210-4900  
Fax: (214) 210-4949  
Web Site: [drew-scientific.com](http://drew-scientific.com)