



# DREW SCIENTIFIC, INC.

## MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: CAT

Product Number: 600064, 600069

Lot: 1664

Expiration Date: 4-15-09

**For use with HEMAVET® and MASCOT™ Hematology Analyzers**

Parameter	Hemavet 800,850,850FS,950,950LV 550,1500,1700,3700, Mascot 850		Hemavet 950 Rev F and above	
	Target	Range	Target	Range
WBC (K/μL)	7.9	7.3 - 8.5	7.9	7.3 - 8.5
NE (K/μL)				
NE/MO/BA (K/μL)				
LY (K/μL)	1.4	0.9 - 1.9	1.4	0.9 - 1.9
MO (K/μL)				--Elev
EO (K/μL)				
BA (K/μL)				--Elev
NRBC (K/μL)				
NE (%)				
NE/MO/BA (%)				
LY (%)	17.4	12.4 - 22.4	17.4	12.4 - 22.4
MO (%)				--Elev
EO (%)				
BA (%)				--Elev
NRBC (%)				
RBC (M/μL)	7.40	7.08 - 7.72	7.40	7.08 - 7.72
HB (g/dL)	13.0	12.4 - 13.6	12.6	12.0 - 13.2
HCT (%)	32.6	28.4 - 37.1	37.6	33.2 - 42.4
MCV (fL)	44.1	40.1 - 48.1	50.9	46.9 - 54.9
MCH (pg)	17.6	16.1 - 19.2	17.1	15.6 - 18.7
MCHC (g/dL)	39.9	33.4 - 48.0	33.6	28.4 - 39.9
RDW (%)	22.6	20.6 - 24.6	22.6	20.6 - 24.6
PLT (K/μL)	191	146 - 236	191	146 - 236

DREW SCIENTIFIC, INC. 4230 Shilling Way Dallas, TX 75237-1023 • USA

### IMPORTANT

Follow mixing instructions from reverse side:  
Wipe threads on bottle and cap before closing

**FOR INVITRO DIAGNOSTIC USE**  
**PROTECT FROM FREEZING**

STORE AT 2° - 8° (35° - 46°)