



DREW SCIENTIFIC, INC.

MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: MOUSE Product Number: 600065, 600070

Lot: 1665-R Expiration Date: 4-15-09

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 1500,1700,3700, Mascot 850		Hemavet 550, Mascot 550	
	Target	Range	Target	Range
WBC (K/μL)	18.9	17.4 - 20.4	18.9	17.4 - 20.4
NE (K/μL)	10.8	8.5 - 13.2		
NE/MO/BA (K/μL)			12.8	10.6 - 15.3
LY (K/μL)	5.8	4.4 - 7.3	5.8	4.4 - 7.3
MO (K/μL)	1.3	0.8 - 1.9		--Elev
EO (K/μL)	0.3	0.0 - 0.9	0.3	0.0 - 0.9
BA (K/μL)	0.8	0.0 - 1.7		--Elev
NRBC (K/μL)				
NE (%)	56.9	48.9 - 64.9		
NE/MO/BA (%)			67.89	60.9 - 74.9
LY (%)	30.6	25.6 - 35.6	30.6	25.6 - 35.6
MO (%)	6.9	4.6 - 9.2		--Elev
EO (%)	1.5	0.0 - 4.5	1.5	0.0 - 4.5
BA (%)	4.1	0.0 - 8.2		--Elev
NRBC (%)				
RBC (M/μL)	8.52	8.13 - 8.91	8.52	8.13 - 8.91
HB (g/dL)	15.2	14.6 - 15.8	15.2	14.6 - 15.8
HCT (%)	45.5	40.2 - 51.1	45.5	40.2 - 51.1
MCV (fL)	53.4	49.4 - 57.4	53.4	49.4 - 57.4
MCH (pg)	17.9	16.4 - 19.5	17.9	16.4 - 19.5
MCHC (g/dL)	33.5	28.6 - 39.4	33.5	28.6 - 39.4
RDW (%)	22.8	20.8 - 24.8	22.8	20.8 - 24.8
PLT (K/μL)	726	566 - 886	726	566 - 886

DREW SCIENTIFIC, INC. 4230 Shilling Way Dallas, TX 75237-1023 • USA

IMPORTANT

Follow mixing instructions from reverse side:
Wipe threads on bottle and cap before closing

FOR INVITRO DIAGNOSTIC USE

PROTECT FROM FREEZING

STORE AT 2° - 8° (35° - 46°)