



# DREW SCIENTIFIC, INC.

## MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: HORSE

Product Number: 600063, 600071

Lot: 1667-R

Expiration Date: 4-15-09

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 1500,1700,3700, Mascot 850		Hemavet 550, Mascot 550	
	Target	Range	Target	Range
WBC (K/μL)	19.1	17.6 - 20.6	19.1	17.6 - 20.6
NE (K/μL)	10.6	8.9 - 12.4		
NE/MO/BA (K/μL)			11.4	9.1 - 13.9
LY (K/μL)	6.8	5.4 - 8.4	6.8	5.4 - 8.4
MO (K/μL)	0.6	0.1 - 1.1		--Elev
EO (K/μL)	0.9	0.0 - 2.0	0.9	0.0 - 2.0
BA (K/μL)	0.3	0.0 - 1.1		--Elev
NRBC (K/μL)				
NE (%)	55.4	50.4 - 60.4		
NE/MO/BA (%)			59.75	51.8 - 67.8
LY (%)	35.6	30.6 - 40.6	35.6	30.6 - 40.6
MO (%)	2.9	0.6 - 5.2		--Elev
EO (%)	4.6	0.0 - 9.6	4.6	0.0 - 9.6
BA (%)	1.5	0.0 - 5.6		--Elev
NRBC (%)				
RBC (M/μL)	8.70	8.31 - 9.09	8.70	8.31 - 9.09
HB (g/dL)	15.8	15.2 - 16.4	15.8	15.2 - 16.4
HCT (%)	39.7	34.6 - 45.1	39.7	34.6 - 45.1
MCV (fL)	45.6	41.6 - 49.6	45.6	41.6 - 49.6
MCH (pg)	18.1	16.7 - 19.7	18.1	16.7 - 19.7
MCHC (g/dL)	39.8	33.7 - 47.4	39.8	33.7 - 47.4
RDW (%)	22.7	20.7 - 24.7	22.7	20.7 - 24.7
PLT (K/μL)	598	438 - 758	598	438 - 758

DREW SCIENTIFIC, INC. 4230 Shilling Way Dallas, TX 75237-1023 • USA

### IMPORTANT

Follow mixing instructions from reverse side:  
Wipe threads on bottle and cap before closing

**FOR INVITRO DIAGNOSTIC USE**

**PROTECT FROM FREEZING**

STORE AT 2° - 8° (35° - 46°)