



# DREW SCIENTIFIC, INC.

## MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: CAT      Product Number: 600064, 600069

Lot: 1674      Expiration Date: 7-15-09

**For use with HEMAVET® and MASCOT™ Hematology Analyzers**

Parameter	Hemavet 800,850,850FS,950,950LV 550,1500,1700,3700, Mascot 850		Hemavet 950, 1700 Rev F and above	
	Target	Range	Target	Range
WBC (K/μL)	8.2	7.6 - 8.8	8.2	7.6 - 8.8
NE (K/μL)				
NE/MO/BA (K/μL)				
LY (K/μL)	1.3	0.8 - 1.8	1.3	0.8 - 1.8
MO (K/μL)				--Elev
EO (K/μL)				
BA (K/μL)				--Elev
NRBC (K/μL)				
NE (%)				
NE/MO/BA (%)				
LY (%)	15.8	10.8 - 20.8	15.8	10.8 - 20.8
MO (%)				--Elev
EO (%)				
BA (%)				--Elev
NRBC (%)				
RBC (M/μL)	7.72	7.40 - 8.04	7.72	7.40 - 8.04
HB (g/dL)	13.1	12.5 - 13.7	13.1	12.5 - 13.7
HCT (%)	36.7	32.2 - 41.4	37.8	33.3 - 42.6
MCV (fL)	47.5	43.5 - 51.5	49.0	45.0 - 53.0
MCH (pg)	17.0	15.5 - 18.5	17.0	15.5 - 18.5
MCHC (g/dL)	35.7	30.2 - 42.5	34.6	29.3 - 41.1
RDW (%)	23.7	21.7 - 25.7	23.7	21.7 - 25.7
PLT (K/μL)	188	143 - 233	188	143 - 233

**DREW SCIENTIFIC, INC. 4230 Shilling Way Dallas, TX 75237-1023 • USA**

### IMPORTANT

Follow mixing instructions from reverse side:  
Wipe threads on bottle and cap before closing

**FOR INVITRO DIAGNOSTIC USE**

**PROTECT FROM FREEZING**

STORE AT 2° - 8° (35° - 46°)