



DREW SCIENTIFIC, INC.

MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: MOUSE

Product Number: 600065, 600070

Lot: 1675

Expiration Date: 7-15-09

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 1500,1700,3700, Mascot 850		Hemavet 550, Mascot 550	
	Target	Range	Target	Range
WBC (K/μL)	19.7	18.2 - 21.2	19.7	18.2 - 21.2
NE (K/μL)	11.2	8.9 - 13.7		
NE/MO/BA (K/μL)			13.3	11.0 - 15.8
LY (K/μL)	6.0	4.6 - 7.5	6.0	4.6 - 7.5
MO (K/μL)	1.4	0.8 - 2.0		--Elev
EO (K/μL)	0.4	0.0 - 1.1	0.4	0.0 - 1.1
BA (K/μL)	0.8	0.0 - 1.7		--Elev
NRBC (K/μL)				
NE (%)	56.8	48.8 - 64.8		
NE/MO/BA (%)			67.70	60.7 - 74.7
LY (%)	30.3	25.3 - 35.3	30.3	25.3 - 35.3
MO (%)	6.9	4.6 - 9.2		--Elev
EO (%)	2.1	0.0 - 5.1	2.1	0.0 - 5.1
BA (%)	4.0	0.0 - 8.1		--Elev
NRBC (%)				
RBC (M/μL)	8.84	8.45 - 9.23	8.84	8.45 - 9.23
HB (g/dL)	15.4	14.8 - 16.0	15.4	14.8 - 16.0
HCT (%)	46.8	41.3 - 52.5	46.8	41.3 - 52.5
MCV (fL)	52.9	48.9 - 56.9	52.9	48.9 - 56.9
MCH (pg)	17.4	16.0 - 18.9	17.4	16.0 - 18.9
MCHC (g/dL)	32.9	28.2 - 38.7	32.9	28.2 - 38.7
RDW (%)	23.9	21.9 - 25.9	23.9	21.9 - 25.9
PLT (K/μL)	701	541 - 861	701	541 - 861

DREW SCIENTIFIC, INC. 4230 Shilling Way Dallas, TX 75237-1023 • USA

IMPORTANT

Follow mixing instructions from reverse side:

Wipe threads on bottle and cap before closing

FOR INVITRO DIAGNOSTIC USE

PROTECT FROM FREEZING

STORE AT 2° - 8° (35° - 46°)