



DREW SCIENTIFIC, INC.

MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: HORSE Product Number: 600063, 600071

Lot: 1677 Expiration Date: 7-15-09

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 1500,1700,3700, Mascot 850		Hemavet 550, Mascot 550	
	Target	Range	Target	Range
WBC (K/μL)	19.4	17.9 - 20.9	19.4	17.9 - 20.9
NE (K/μL)	10.8	9.1 - 12.7		
NE/MO/BA (K/μL)			11.5	9.2 - 14.1
LY (K/μL)	6.8	5.4 - 8.4	6.8	5.4 - 8.4
MO (K/μL)	0.5	0.1 - 1.0		--Elev
EO (K/μL)	1.0	0.1 - 2.2	1.0	0.1 - 2.2
BA (K/μL)	0.2	0.0 - 1.0		--Elev
NRBC (K/μL)				
NE (%)	55.9	50.9 - 60.9		
NE/MO/BA (%)			59.50	51.5 - 67.5
LY (%)	35.2	30.2 - 40.2	35.2	30.2 - 40.2
MO (%)	2.7	0.4 - 5.0		--Elev
EO (%)	5.3	0.3 - 10.3	5.3	0.3 - 10.3
BA (%)	0.9	0.0 - 5.0		--Elev
NRBC (%)				
RBC (M/μL)	9.14	8.75 - 9.53	9.14	8.75 - 9.53
HB (g/dL)	15.9	15.3 - 16.5	15.9	15.3 - 16.5
HCT (%)	41.6	36.3 - 47.2	41.6	36.3 - 47.2
MCV (fL)	45.5	41.5 - 49.5	45.5	41.5 - 49.5
MCH (pg)	17.4	16.1 - 18.9	17.4	16.1 - 18.9
MCHC (g/dL)	38.2	32.4 - 45.4	38.2	32.4 - 45.4
RDW (%)	23.8	21.8 - 25.8	23.8	21.8 - 25.8
PLT (K/μL)	624	464 - 784	624	464 - 784

DREW SCIENTIFIC, INC. 4230 Shilling Way Dallas, TX 75237-1023 • USA

IMPORTANT

Follow mixing instructions from reverse side:
Wipe threads on bottle and cap before closing

FOR INVITRO DIAGNOSTIC USE

PROTECT FROM FREEZING

STORE AT 2^o - 8^o (35^o - 46^o)