



DREW SCIENTIFIC, INC.

MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: CAT Product Number: 600064, 600069

Lot 1694 Expiration Date: 01-15-2010

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 550,1500,1700,3700, Mascot 850		Hemavet 950, 1700 Rev F and above	
	Target	Range	Target	Range
WBC (K/μL)	8.0	7.4 - 8.6	8.0	7.4 - 8.6
NE (K/μL)				
NE/MO/BA (K/μL)				
LY (K/μL)	1.5	1.0 - 2.0	1.5	1.0 - 2.0
MO (K/μL)				--Elev
EO (K/μL)				
BA (K/μL)				--Elev
NRBC (K/μL)				
NE (%)				
NE/MO/BA (%)				
LY (%)	18.4	13.4 - 23.4	18.4	13.4 - 23.4
MO (%)				--Elev
EO (%)				
BA (%)				--Elev
NRBC (%)				
RBC (M/μL)	7.33	7.01 - 7.65	7.33	7.01 - 7.65
HB (g/dL)	12.3	11.7 - 12.9	12.3	11.7 - 12.9
HCT (%)	32.8	28.6 - 37.3	36.2	31.8 - 40.8
MCV (fL)	44.8	40.8 - 48.8	49.4	45.4 - 53.4
MCH (pg)	16.8	15.3 - 18.4	16.8	15.3 - 18.4
MCHC (g/dL)	37.4	31.3 - 45.1	34.0	28.6 - 40.5
RDW (%)	23.0	21.0 - 25.0	23.0	21.0 - 25.0
PLT (K/μL)	189	144 - 234	189	144 - 234

DREW SCIENTIFIC, INC. 4230 Shilling Way Dallas, TX 75237-1023 • USA

IMPORTANT

Follow mixing instructions from reverse side:
Wipe threads on bottle and cap before closing

FOR INVITRO DIAGNOSTIC USE

PROTECT FROM FREEZING

STORE AT 2° - 8° (35° - 46°)