



DREW SCIENTIFIC, INC.

MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: DOG Product Number: 600066, 600068

Lot 1696 Expiration Date: 01-15-2010

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 1500,1700,3700, Mascot 850		Hemavet 550, Mascot 550	
	Target	Range	Target	Range
WBC (K/μL)	12.3	11.7 - 12.9	12.3	11.7 - 12.9
NE (K/μL)	9.2	7.8 - 10.6		
NE/MO/BA (K/μL)			9.5	8.1 - 11.0
LY (K/μL)	1.4	0.8 - 2.1	1.4	0.8 - 2.1
MO (K/μL)	0.1	0.0 - 0.4		--Elev
EO (K/μL)	1.4	0.8 - 2.2	1.4	0.8 - 2.2
BA (K/μL)	0.2	0.0 - 0.5		--Elev
NRBC (K/μL)				
NE (%)	74.3	66.3 - 82.3		
NE/MO/BA (%)			76.80	68.8 - 84.8
LY (%)	11.5	6.5 - 16.5	11.5	6.5 - 16.5
MO (%)	0.8	0.0 - 3.1		--Elev
EO (%)	11.7	6.7 - 16.7	11.7	6.7 - 16.7
BA (%)	1.7	0.0 - 3.7		--Elev
NRBC (%)				
RBC (M/μL)	5.90	5.61 - 6.19	5.90	5.61 - 6.19
HB (g/dL)	14.8	14.2 - 15.4	14.8	14.2 - 15.4
HCT (%)	37.9	33.8 - 42.2	37.9	33.8 - 42.2
MCV (fL)	64.2	60.2 - 68.2	64.2	60.2 - 68.2
MCH (pg)	25.0	22.8 - 27.3	25.0	22.8 - 27.3
MCHC (g/dL)	38.9	33.5 - 45.4	38.9	33.5 - 45.4
RDW (%)	18.7	16.7 - 20.7	18.7	16.7 - 20.7
PLT (K/μL)	433	368 - 498	433	368 - 498

DREW SCIENTIFIC, INC. 4230 Shilling Way Dallas, TX 75237-1023 • USA

IMPORTANT

Follow mixing instructions from reverse side:
Wipe threads on bottle and cap before closing

FOR INVITRO DIAGNOSTIC USE

PROTECT FROM FREEZING

STORE AT 2° - 8° (35° - 46°)