



# DREW SCIENTIFIC, INC.

## MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: HORSE Product Number: 600063, 600071

Lot 1697 Expiration Date: 01-15-2010

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 1500,1700,3700, Mascot 850		Hemavet 550, Mascot 550	
	Target	Range	Target	Range
WBC (K/μL)	19.2	17.7 - 20.7	19.2	17.7 - 20.7
NE (K/μL)	6.7	5.3 - 8.3		
NE/MO/BA (K/μL)			7.4	5.4 - 9.6
LY (K/μL)	7.7	6.2 - 9.3	7.7	6.2 - 9.3
MO (K/μL)	0.1	0.0 - 0.6		--Elev
EO (K/μL)	4.2	2.9 - 5.5	4.2	2.9 - 5.5
BA (K/μL)	0.5	0.0 - 1.4		--Elev
NRBC (K/μL)				
NE (%)	34.9	29.9 - 39.9		
NE/MO/BA (%)			38.41	30.4 - 46.4
LY (%)	40.0	35.0 - 45.0	40.0	35.0 - 45.0
MO (%)	0.7	0.0 - 3.0		--Elev
EO (%)	21.6	16.6 - 26.6	21.6	16.6 - 26.6
BA (%)	2.9	0.0 - 7.0		--Elev
NRBC (%)				
RBC (M/μL)	8.73	8.34 - 9.12	8.73	8.34 - 9.12
HB (g/dL)	15.0	14.4 - 15.6	15.0	14.4 - 15.6
HCT (%)	38.4	33.4 - 43.8	38.4	33.4 - 43.8
MCV (fL)	44.0	40.0 - 48.0	44.0	40.0 - 48.0
MCH (pg)	17.2	15.8 - 18.7	17.2	15.8 - 18.7
MCHC (g/dL)	39.0	32.9 - 46.7	39.0	32.9 - 46.7
RDW (%)	23.9	21.9 - 25.9	23.9	21.9 - 25.9
PLT (K/μL)	596	436 - 756	596	436 - 756

DREW SCIENTIFIC, INC. 4230 Shilling Way Dallas, TX 75237-1023 • USA

### IMPORTANT

Follow mixing instructions from reverse side:  
Wipe threads on bottle and cap before closing

**FOR INVITRO DIAGNOSTIC USE**

**PROTECT FROM FREEZING**

STORE AT 2° - 8° (35° - 46°)